

CAVITY/VDRL/PRINTED/ADHESIVE SLIDES

**BLUE
STAR**

PRICE LIST
APRIL 2026

BLUE STAR

1. BLUE STAR (BLOOD SMEAR SLIDES)

Size: 75 mm x 25 mm x 1.35 mm thick, Edges Rounded

Per Pkt. of 10 Nos.

Rs.160/-

2. BLUE STAR (MICRO CONCAVITY SLIDES)

Size: 75 mm x 25 mm x 1.35mm thick

Polished spherical concavities are 15 mm in diameter x 0.5 mm deep with beveled and polished edges

(a) One Concavity

(b) Two Concavities

(c) Three Concavities

Rs.150/-

Rs.155/-

Rs.160/-

3. BLUE STAR V.D.R.L. SLIDE (KLINE CONCAVITY SLIDES)

each Size: 75 mm x 56 mm x 3 mm thick with spherical concavities, each 16 mm in diameter x 1.75 mm deep with beveled and polished edges 2, 4, 6, 8, 10 & 12 Concavities

Per each

Rs.150/-

4. BLUE STAR BLOOD GROUPING SLIDES.

Size: 100 mm x 100 mm x 10 mm thick with 12 Concavities,

Per each 21 mm in diameter x 4 mm deep. Surface of the slide is frosted to provide smooth writing. Surface edges and corners are beveled and polished

Rs.600/-

5. BLUE STAR PRINTED SLIDES

Size: 75 mm X 25mm X 1.10 mm, One end has 19mm.

printed surface on one side for marking in White Coating

Per Pkt of 20 Nos.

Rs.140/-

6. BLUE STAR ADHESIVE SLIDES.

Size: 75 mm X 25mm X 1.10 mm, One end has 19mm

Special coating which enhance adhesion

(a) Positive Charge Slides


(b) Polysine Slides

(c) Silane Slides

(d) Hydrophile Slides

(e) Cytology Slide

Rs.450/-



POLAR INDUSTRIAL CORPORATION

Shree Laxmi Woollen Mills Compound, Off Dr. E. Moses Road,

Mahalaxmi, Mumbai - 400 011 (India)

Phone: +91 2249731752 / 9324148397

Email: polar1bluestar@gmail.com

Website: <http://www.bluestarslides.com>

**BLUE
STAR**



MICRO SLIDES & COVER GLASSES



BLUE STAR

MICRO SLIDES

Manufactured from best, selected and optically flat Glass Sheets. Micro slides are uniform in size with ground edges most suitable for Microscopic work.

Size: 75mm X 25mm	Per Pkt of 50 Nos.
1. BLUE STAR (PIC 2) Thickness: 1.45 mm	Rs.79/-
2. BLUE STAR (PIC 1) Thickness: 1.35 mm	Rs.80/-
3. BLUE STAR (PIC 1 PLUS) Thickness: 1.35 mm (Interleaved)	Rs.85/-
4. BLUE STAR (DELUXE) Thickness: 1.25 mm	Rs.90/-
5. BLUE STAR (SUPER DELUXE) Thickness: 1.10 mm	Rs.100/-
6. BLUE STAR (SPECIAL SUPER DELUXE) Thickness: 1.00 mm	Rs.105/-
7. BLUE STAR (FROSTED MICRO SLIDES)	
Thickness: 1.35 mm	Rs.105/-
Thickness: 1.25 mm	Rs.110/-
Thickness: 1.10 mm	Rs.120/-
Thickness: 1.00 mm	Rs.125/-
One end has 19 mm fine frosted surface on both the sides for marking	
8. BLUE STAR (SPECIAL SIZE)	
Size: 75 mm x 50 mm	
Thickness: 1.35 mm	Rs.270/-
Thickness: 1.10 mm	Rs.320/-

PACKING DETAILS:

Our Standard Packing is 50 Nos. in a Packet duly gelatined, then 10 Packets in an inner corrugated box sealed with BLUE STAR Tape & likewise 10 corrugated boxes in a strong printed outer corrugated box sealed with Blue & White BLUE STAR Tape. i.e. Packing of 100 Packets in a corrugated box.

- N. B. 1) Odd sizes for all the items will be quoted and supplied on request.
 2) Our all products are pre cleaned and ready to use.
 3) Variation in Size + 1mm and Thickness + 0.1mm.
 4) Prices shown are maximum. Local Taxes extra.



BLUE STAR

COVER GLASSES

Manufactured from IMPORTED mechanically drawn glass, which is superior than mouth blown glass. Cover glasses are optically flat & free from waviness.

No.0 0.08 mm to 0.13 mm 1/250"
 No.1 0.13 mm to 0.16 mm 1/200"

Available Sizes	Per Pkt. of 10 grams
	No.1 No.0
1. SQUARE 18 mm x 18 mm, 19 mm x 19 mm, 22 mm x 22 mm, 24 mm x 24 mm, 25 mm x 25 mm	Rs.105/- Rs.130/-
2. RECTANGULAR 22 mm x 25 mm, 22 mm x 30 mm, 22 mm x 40 mm, 22 mm x 50 mm	Rs.107/- Rs.135/-
3. RECTANGULAR (SPECIAL) 22 mm x 60 mm, 24 mm x 30 mm, 24 mm x 40 mm, 24 mm x 50 mm, 24 mm x 60 mm, 25 mm x 50 mm	Rs.113/- Rs.140/-
4. CIRCULAR 18 mm Circular, 19 mm Circular, 22 mm Circular, 24 mm Circular	Rs.170/- Rs.200/-
5. CIRCULAR (SPECIAL) 12 mm Circular, 16 mm Circular	Rs.180/- Rs.210/-
6. HAEMOCYTOMETER COVER GLASS 20 mm x 25 mm (No.4 i.e. 0.4mm Thick) Packing 20 Nos. in a packet	Rs.300/-

PACKING DETAILS:

Our Standard Packing is 10 grams in Ivory colour attractive Plastic Packet and likewise 20 packets in a very attractive laminated folding box. Packed dry & perfectly gelatined to avoid moisture.

CAUTION NOTICE

This is to draw to your attention that some Persons/Dealers are packing & selling spurious & inferior quality Micro slides & Cover glasses in packing exactly similar to ours bearing our Regd. Trade Mark "BLUE STAR". These are not genuine. In view of this we suggest to our valued patrons to procure "BLUE STAR" Micro slides & Cover glasses either directly from us or from our Authorised Stockists/Distributors.