



**Laboratory  
Glassware**  
Mandatory in India

## PRICE LIST 2025-26



**MSME**  
सूक्ष्म, लघु एवं मध्यम उद्यम  
MICRO, SMALL & MEDIUM ENTERPRISES



एन एस आई सी  
**NSIC**



**An ISO 9001:2015 Certified Company**

**Manufactures & Suppliers of : All type of Laboratory Glassware & Lab Equipments**



[www.ashaglassware.com](http://www.ashaglassware.com)



[ashascientificworks@ymail.com](mailto:ashascientificworks@ymail.com)

# COMPANY PROFILE & CERTIFICATIONS

**An ISO 9001:2015 Certified Company**

**Manufactures & Suppliers of : All type of Laboratory Glassware & Lab Equipments**



[www.ashaglassware.com](http://www.ashaglassware.com)



[ashasientificworks@ymail.com](mailto:ashasientificworks@ymail.com)



# COMPANY PROFILE & CERTIFICATIONS

**Asha Scientific Works** based at Ambala (INDIA) is one of the trusted and fastest growing Manufacturer & Exporters of best quality Scientific Laboratory Glassware & Instruments in the country. Products manufactured by us are one of the finest qualities to add durability & strength. To meet the different demands of our valuable customers, we offer Laboratory Glassware & Instruments in various sizes, types and specifications. We specialize in manufacturing Glassware as per customer's requirements and drawings. With the team of best professionals and innovative technologies, we ensure the right time delivery of our products. Being a major player in the industry, the zeal still prevails in our employees towards customer satisfaction. We empower clients by offering optimum quality Laboratory Equipment within the stipulated time frame at affordable prices. Our products are extensively used by the Chemical Industries, Laboratories, Pharmaceutical Industry, Petrochemicals, Refineries, Specialty Chemicals, Healthcare, Educational Segments, Institutions, Research & QC Laboratories, Agrochemical, Renewable Energy, Food Safety, Medical and Environmental Monitoring.

We understand that your lab Equipment / Instruments plays a crucial role in your work. That's why we offer more than just products; we provide solutions. If you are looking for high quality Scientific Laboratory Glassware, Equipment & Instruments to support your lab's research and experiments, look no further. Contact **Asha Scientific Works** today and let us assist you in finding the perfect solutions for your lab needs.





# FLASKS

## A19

### Boiling Flasks, Round Bottom, ISI Marked Single Neck

- Made from Borosilicate 3.3 Heat Resistant Glass.
- Compliance with IS 1381 Part -1

IS : 1381 Part-1



CM/L-9600102615



Cat. No.	Capacity (mL)	Pack of
A19/1	25	10
A19/2	50	6
A19/3	100	6
A19/4	150	6
A19/5	250	6
A19/6	500	4
A19/7	1000	2
A19/8	2000	2
A19/9	3000	2
A19/10	5000	2
A19/11	10000	1

## A20

### Boiling Flasks, Flat Bottom, ISI Marked, Single Neck

- Made from Borosilicate 3.3 Heat Resistant Glass.
- Compliance with IS 1381 Part -1

IS : 1381 Part-1



CM/L-9600102615

Cat. No.	Capacity (mL)	Pack of
A20/1	25	10
A20/2	50	6
A20/3	100	6
A20/4	150	6
A20/5	250	6
A20/6	500	4
A20/7	1000	2
A20/8	2000	2
A20/9	3000	2
A20/10	5000	2
A20/11	10000	1





## BEAKERS

IS : 2619



CM/L -9600102514

### A11 ISI Marked Beakers Low Form with Spout

- Made from Borosilicate 3.3 Heat Resistant Glass.

Cat. No.	Capacity (mL)	Diameter (mm)	Height (mm)	Pack of
A11/1	5	22	30	10
A11/2	10	26	35	10
A11/3	25	34	51	10
A11/4	50	42	58	6
A11/5	100	50	70	6
A11/6	150	60	88	6
A11/7	250	68	95	6
A11/8	400	80	110	6
A11/9	500	83	118	6
A11/10	600	90	125	6
A11/11	800	100	140	6
A11/12	1000	105	157	6
A11/13	2000	132	185	4
A11/14	3000	152	210	1
A11/15	5000	170	270	1
A11/16	10000*	220	350	1

\* Marked is Non-ISI Product.

### A16 Flasks, Erlenmeyer, ISI Marked , Narrow Mouth

- Made from Borosilicate 3.3 Heat Resistant Glass.
- Compliance with IS 1381 Part-1

IS : 1381 Part-1



CM/L-9600102615



## FLASKS

Cat. No.	Capacity (mL)	O.D. × Height (approx. mm)	Neck O.D. (approx. mm)	Pack of
A16/1	25	42 × 75	22	1
A16/2	50	51 × 90	22	6
A16/3	100	64 × 105	26	6
A16/4	150	74 × 118	28	6
A16/5	250	85 × 145	34	6
A16/6	500	105 × 180	34	6
A16/7	1000	131 × 220	42	2
A16/8	2000	166 × 280	48	1
A16/9	3000	187 × 310	48	1
A16/10	5000	220 × 365	50	1







# VOLUMETRIC FLASKS

## A28

### Volumetric Flasks, ISI Marked, Clear, Class 'A'

- Compliance with IS : 915 standard.
- With individual Class 'A' Certificate.
- Printed in Anti Acid White/Amber/Blue Colour.
- Come with PP/Glass Penny Stopper.
- Calibrated at 27°C.

IS : 915



CM/L-9600125013

Cat. No.	Capacity (mL)	Neck Size	Tolerance (±mL)	Pack of
A28/1	5	7/16	0.02	1
A28/2	10	7/16	0.025	1
A28/3	20	10/19	0.040	1
A28/	25	10/19	0.040	1
A28/	50	12/21	0.060	1
A28/	100	12/21	0.10	1
A28/7	200	14/23	0.150	1
A28/	250	14/23	0.150	1
A28/9	500	19/26	0.25	1
A28/1	1000	24/29	0.400	1
A28/1	2000	29/32	0.600	1
A28/1	5000	34/35	1.200	1



IS : 915



CM/L-9600125013

## A29

### Volumetric Flasks, ISI Marked, Clear, Class 'B'

- Compliance with IS : 915 standard.
- Accuracy as per Class 'B'.
- Printed in Anti Acid White/Amber/Blue Colour.
- Come with PP/Glass Penny Stopper.
- Calibrated at 27°C.

Cat. No.	Capacity (mL)	Neck Size	Tolerance (±mL)	Pack of
A29/	5	7/16	0.05	2
A29/	10	7/16	0.05	2
A29/	20	10/19	0.08	2
A29/	25	10/19	0.08	2
A29/	50	12/21	0.12	2
A29/	100	14/23	0.20	2
A29/7	200	14/23	0.30	2
A29/8	250	14/23	0.30	2
A29/9	500	19/26	0.50	2
A29/1	1000	24/29	0.80	2
A29/1	2000	29/32	1.20	1
A29/1	5000	34/35	2.50	1





# VOLUMETRIC FLASKS

## A30

### Volumetric Flasks, ISI Marked, Amber, Class 'A'

- Compliance with IS : 915 standard.
- With individual Class 'A' Certificate.
- Printed in Anti Acid White Colour.
- Come with PP/Glass Stopper.
- Calibrated at 27°C.

IS : 915



CM/L-9600125013



Cat. No.	Capacity (mL)	Neck Size	Tolerance ( $\pm$ mL)	Pack of
A30/1	5	7/16	0.025	1
A30/2	10	7/16	0.025	1
A30/3	20	10/19	0.040	1
A30/4	25	10/19	0.040	1
A30/5	50	12/21	0.060	1
A30/6	100	14/23	0.100	1
A30/7	200	14/23	0.150	1
A30/8	250	14/23	0.150	1
A30/9	500	19/26	0.250	1
A30/10	1000	24/29	0.400	1
A30/11	2000	29/32	0.600	1
A30/12	5000	34/35	1.200	1

## A31

### Volumetric Flasks, ISI Marked, Amber, Class 'B'

- Compliance with IS : 915 standard.
- Accuracy as per Class 'B'.
- Printed in Anti Acid White Colour.
- Come with PP/Glass Stopper.
- Calibrated at 27°C.

IS : 915




CM/L-9600125013

Cat. No.	Capacity (mL)	Neck Size	Tolerance ( $\pm$ mL)	Pack of
A31/1	5	10/19	0.05	2
A31/2	10	10/19	0.05	2
A31/3	20	10/19	0.08	2
A31/4	25	10/19	0.08	2
A31/5	50	12/21	0.12	2
A31/6	100	14/23	0.20	2
A31/7	200	14/23	0.30	2
A31/8	250	14/23	0.30	2
A31/9	500	19/26	0.50	2
A31/10	1000	24/29	0.80	2
A31/11	2000	29/32	1.20	1
A31/12	5000	34/35	2.50	1





# CYLINDERS



## A42

### Measuring Cylinders **ISI Marked**, Hexagonal Base,

- Compliance with IS : 878 Standard.
- With individual Class 'A' Certificate.
- Printed in Anti Acid White Colour.
- Highly resistant to chemical attack.
- Calibrated at 27°C (TC).

**Class 'A'**

IS : 878



CM/L-9600115717

Cat. No.	Capacity (mL)	Graduation (mL)	Tolerance (±mL)	Pack of
A42/1	5	0.1	0.05	1
A42/2	10	0.2	0.10	1
A42/3	25	0.5	0.25	1
A42/4	50	1.0	0.50	1
A42/5	100	1.0	0.50	1
A42/6	250	2.0	1.00	1
A42/7	500	5.0	2.50	1
A42/8	1000	10.0	5.00	1
A42/9	2000	20.0	10.00	1

## A43

### Measuring Cylinders **ISI Marked**, Hexagonal Base,

- Compliance with IS : 878 Standard.
- Accuracy as per Class 'B'.
- Printed in Anti Acid White Colour.
- Highly resistant to chemical attack.
- Calibrated at 27°C (TC).

**Class 'B'**

IS : 878



CM/L-9600115717


Cat. No.	Capacity (mL)	Graduation (mL)	Tolerance (±mL)	Pack of
A43/1	5	0.1	0.10	2
A43/2	10	0.2	0.20	2
A43/3	25	0.5	0.50	2
A43/4	50	1.0	1.00	2
A43/5	100	1.0	1.00	2
A43/6	250	2.0	2.00	1
A43/7	500	5.0	5.00	1
A43/8	1000	10.0	10.00	1
A43/9	2000	20.0	20.00	1







# CYLINDERS



## A40

### Measuring Cylinders, **ISI Marked**, Round Base, Class 'A'

- Compliance with IS : 878 Standard.
- With individual Class 'A' Certificate.
- Printed in Anti Acid White Colour.
- Highly resistant to chemical attack.
- Calibrated at 27°C (TC).

IS : 878



CM/L-9600115717

Cat. No.	Capacity (mL)	Graduation (mL)	Tolerance ( $\pm$ mL)	Pack of
A40/1	5	0.1	0.05	1
A40/2	10	0.2	0.10	1
A40/3	25	0.5	0.25	1
A40/4	50	1.0	0.50	1
A40/5	100	1.0	0.50	1
A40/6	250	2.0	1.00	1
A40/7	500	5.0	2.50	1
A40/8	1000	10.0	5.00	1
A40/9	2000	20.0	10.00	1

## A41

### Measuring Cylinders, **ISI Marked**, Round Base, Class 'B'


- Compliance with IS : 878 Standard.
- Accuracy as per Class 'B'.
- Printed in Anti Acid White Colour.
- Highly resistant to chemical attack.
- Calibrated at 27°C (TC).

IS : 878



CM/L-9600115717

Cat. No.	Capacity (mL)	Graduation (mL)	Tolerance ( $\pm$ mL)	Pack of
A41/1	5	0.1	0.10	2
A41/2	10	0.2	0.20	2
A41/3	25	0.5	0.50	2
A41/4	50	1.0	1.00	2
A41/5	100	1.0	1.00	2
A41/6	250	2.0	2.00	2
A41/7	500	5.0	5.00	1
A41/8	1000	10.0	10.00	1
A41/9	2000	20.0	20.00	1





# PIPETTES

## A51

### Volumetric Pipettes, ISI Marked One Mark, Class 'A'

- Compliance with IS : 1117 standard.
- With individual Class 'A' Certificate.
- Each pipette with unique identification number.
- Calibrated and certified at 27°C.
- Calibrated for delivery (TD, Ex).

IS : 1117



Cat. No.	Capacity (mL)	Tolerance ( $\pm$ mL)	Pack of
A51/1	1	0.006	10
A51/2	2	0.010	10
A51/3	5	0.015	10
A51/4	10	0.020	10
A51/5	20	0.030	10
A51/6	25	0.030	10
A51/7	50	0.050	10
A51/8	100	0.080	10

Pipettes with Safety Bulb cost 10% extra over list price.



## A52

### Volumetric Pipettes ISI Marked, One Mark, Class 'B'

- Compliance with IS : 1117 standard.
- Accuracy as per Class 'B'.
- Calibrated and certified at 27°C.
- Calibrated for delivery (TD, Ex).

IS : 1117



Cat. No.	Capacity (mL)	Tolerance ( $\pm$ mL)	Pack of
A52/1	1	0.015	10
A52/2	2	0.020	10
A52/3	5	0.030	10
A52/4	10	0.040	10
A52/5	20	0.060	10
A52/6	25	0.060	10
A52/7	50	0.100	10
A52/8	100	0.150	10

Pipettes with Safety Bulb cost 10% extra over list price.





**Laboratory  
Glassware**  
**Mandatory in India**

**“Quality Above  
All in Science”**



**Workshop & Office : 2 Lok Nayak Society, Near Housing Board Colony, Ambala Cantt-133001 Haryana, INDIA**  
**Corporate Office : Plot No. 121, Vikas Puri, Sai Baba Mandir, Sarsehri, Ambala-133001 Haryana, INDIA**  
**Phones: +91-98961-51568, +91-74090-02244, +91-98963-80089**



[www.ashaglassware.com](http://www.ashaglassware.com)



[ashascientificworks@ymail.com](mailto:ashascientificworks@ymail.com)