

1. BLUE STAR (BLOOD SMEAR SLIDES)

Size: 75 mm x 25 mm x 1.35 mm thick, Edges Rounded

Per Pkt. of 10 Nos.

Rs.145/-

2. BLUE STAR (MICRO CONCAVITY SLIDES)

Size: 75 mm x 25 mm x 1.45mm thick

Polished spherical concavities are 15 mm in diameter
x 0.5 mm deep with beveled and polished edges

- (a) One Concavity
- (b) Two Concavities
- (c) Three Concavities

Rs.135/-

Rs.140/-

Rs.145/-

3. BLUE STAR V.D.R.L. SLIDE (KLINE CONCAVITY SLIDES)Size: 75 mm x 56 mm x 3 mm thick with spherical concavities, each
16 mm in diameter x 1.75 mm deep with beveled and polished edges
2, 4, 6, 8, 10 & 12 Concavities

Per Each

Rs.135/-

4. BLUE STAR BLOOD GROUPING SLIDES.Size: 100 mm x 100 mm x 10 mm thick with 12 Concavities,
each 21 mm in diameter x 4 mm deep. Surface of the slide
is frosted to provide smooth writing. Surface edges and
corners are beveled and polished

Rs.550/-

5. BLUE STAR PRINTED SLIDESSize: 75 mm X 25mm X 1.10 mm, One end has 19mm
printed surface on one side for marking in White Coating

Per Pkt of 20 Nos.

Rs.125/-

6. BLUE STAR ADHESIVE SLIDESSize: 75 mm X 25mm X 1.10 mm, One end has 19mm
Special coating which enhance adhesion

- (a) Positive Charge Slides
- (b) Polysine Slides
- (c) Silane Slides
- (d) Hydrophile Slides
- (e) Cytology Slide

Rs.400/-

**POLAR INDUSTRIAL CORPORATION**

Shree Laxmi Woollen Mills Compound, Off Dr. E. Moses Road,
Mahalaxmi, Mumbai - 400 011 (India)
Phone: +91 24921752/49731752
Email: polar1bluestar@gmail.com
Website: <http://www.bluestarslides.com>



BLUE STAR®
POLAR INDUSTRIAL CORPORATION





BLUE STAR MICRO SLIDES

Manufactured from best, selected and optically flat Glass Sheets. Micro slides are uniform in size with ground edges most suitable for Microscopic work.

Size: 75mm X 25mm	Per Pkt of 50 Nos.
1. BLUE STAR (PIC 2) Thickness: 1.35 mm	Rs.74/-
2. BLUE STAR (PIC 1) Thickness: 1.35 mm	RS.75/-
3. BLUE STAR (PIC 1 PLUS) Thickness: 1.35 mm (Interleaved)	Rs.80/-
4. BLUE STAR (DELUXE) Thickness: 1.25 mm	Rs.85/-
5. BLUE STAR (SUPER DELUXE) Thickness: 1.10 mm	Rs.95/-
6. BLUE STAR (SPECIAL SUPER DELUXE) Thickness: 1.00 mm	Rs.100/-
7. BLUE STAR (FROSTED MICRO SLIDES) Thickness: 1.35 mm Thickness: 1.25 mm Thickness: 1.10 mm Thickness: 1.00 mm One end has 19 mm fine frosted surface on both the sides for marking	Rs.100/- Rs.105/- Rs.115/- Rs.120/-
8. BLUE STAR (SPECIAL SIZE) Size: 75 mm x 50 mm Thickness: 1.35 mm Thickness: 1.10 mm	Rs.250/- Rs.300/-

PACKING DETAILS:

Our Standard Packing is 50 Nos. in a Packet duly gelatined, then 10 Packets in an inner corrugated box sealed with BLUE STAR Tape & likewise 10 corrugated boxes in a strong printed outer corrugated box sealed with Blue & White BLUE STAR Tape. i.e. Packing of 100 Packets in a corrugated box.

N. B. 1) Odd sizes for all the items will be quoted and supplied on request.

2) Our all products are pre cleaned and ready to use.

3) Variation in Size ± 1 mm and Thickness ± 0.1 mm.

4) Prices shown are maximum. Local Taxes extra.

BLUE STAR COVER GLASSES



Manufactured from IMPORTED mechanically drawn glass, which is superior than mouth blown glass. Cover glasses are optically flat and free from waviness.

No. 0 0.08 mm to 0.13 mm 1/250"
No. 1 0.13 mm to 0.16 mm 1/200"

Available Sizes

1. SQUARE

18 mm x 18 mm, 19 mm x 19 mm,
22 mm x 22 mm, 24 mm x 24 mm,
25 mm x 25 mm

2. RECTANGULAR

22 mm x 25 mm, 22 mm x 30 mm,
22 mm x 40 mm, 22 mm x 50 mm,

3. RECTANGULAR (SPECIAL)

22 mm x 60 mm, 24 mm x 40 mm,
24 mm x 50 mm, 24 mm x 60 mm,
25 mm x 50 mm

4. CIRCULAR

18 mm Circular, 19 mm Circular,
22 mm Circular, 24 mm Circular

5. CIRCULAR (SPECIAL)

12 mm Circular, 16 mm Circular

6. HAEMOCYTOMETER COVER GLASS

20 mm x 25 mm (No.4 i.e. 0.4mm Thick)
Packing 20 Nos. in a packet

Per Pkt. of 10 grams

No.1 No.0
Rs.95/- Rs.120/-

Rs.97/- Rs.125/-

Rs.103/- Rs.130/-

Rs.155/- Rs.185/-

Rs.165/- Rs.195/-

Rs.270/-

PACKING DETAILS:

Our Standard Packing is 10 grams in Ivory colour attractive Plastic Packet and likewise 20 Packets in a very attractive laminated folding box. Packed dry & perfectly gelatined to avoid moisture.

CAUTION NOTICE

This is to draw to your attention that some Persons/Dealers are packing & selling spurious & inferior quality Micro slides & Cover glasses in packing exactly similar to ours bearing our Regd. Trade Mark "BLUE STAR". These are not genuine. In view of this we suggest to our valued patrons to procure "BLUE STAR" Micro slides & Cover glasses either directly from us or from our Authorised Stockists/Distributors.